



SHOP FROM HOME Hearing Aid Catalog

FDA LISTED
100% DIGITAL
AMERICAN-MADE TECHNOLOGY



Hearing aid technology meticulously engineered
to deliver **better sound** and **improved hearing!**

 **HearingHelp express**
an IntriCon company

TRUSTED BY 750,000+ CUSTOMERS SINCE 1979

INDEX



About Hearing Help Express	3
About our hearing aids	4
Our process	5
Choosing the right hearing aid	6
Explanation of product features	7
EarMate® Series hearing aids	8-10
EarMate-4220	
EarMate-4330	
EarMate-4210	
Apollo® Series hearing aids	11-13
Apollo-6205	
Apollo-E600	
Apollo-6670	
Felix™ Accu-Fit	14-15
Felix-8750	
Ranger® Series hearing aids	16-18
Ranger-8720	
Ranger-7800	
Ranger-7910	
Service plans	19
Hearing aid batteries	20
Helpful accessories	21-24
Order form	25-26
Commonly asked questions	27
Hearing industry survey results	27



TRUSTED BY 750,000+
CUSTOMERS SINCE 1979

1-800-221-2099



Welcome to

HearingHelp express

an IntriCon company

- > High-quality, affordable hearing aids
- > Shop safely from home
- > 45-day risk-free home trial
- > Free support from licensed professionals
- > Payment plans available
- > FREE shipping

About Us

Hearing Help Express's journey began in 1979 to improve the quality of life of those suffering from hearing loss. Hearing Help Express is the oldest direct to consumer hearing aid provider with more than 750,000 satisfied customers. Better yet, with our on-demand system and teleaudiology, you'll get a hearing assessment and recommendation by a trained hearing aid professional as well as a lifetime of hearing help support, all from the comfort of your own home.

Hearing Help Express is your best option for affordable, U.S.-Made hearing aid technology, and you save thousands without sacrificing hearing improvement.

We're ready to help you hear better.



Quality

YOU CAN TRUST & TRY



On-Demand SUPPORT



Teleaudiology

REMOTE CARE

About Our Hearing Aids

Hearing Help Express is owned by Intricon, a U.S. manufacturer of FDA-listed hearing aids. This allows us to offer you the same high-quality hearing aids sold by audiologists and ENT's, keeping quality high and prices low, passing the savings on to you.



- > 100% digital
- > FDA-listed devices
- > American-Made technology

Hearing aid technology meticulously engineered to deliver **better sound** and **improved hearing!**



PROVEN
QUALITY



ACCURATE HEARING
RESULTS



BEST VALUE
HEARING

"I have tried a number of hearing aids including the \$5,000 ones and yours are the best I have ever used!"

Margaret E. – AL

"I have been wearing hearing aids for 35 years and these are the best I have ever had."

William J – MI

"You have great products. After wearing hearing aids for 36 years of all kinds and costs, yours are the best ever."

Col. Paul G - CA

The Hearing Help Express process and how it can work for you

Choose from 4 easy option to help you determine your level of hearing loss



Complete the free
hearing assessment
enclosed

Call 1-800-221-2099

Have your hearing assessed over the phone by a licensed hearing professional

Take our free online hearing assessment
hearinghelp.com/hearing-check

Send us your hearing test (audiogram) by email: info@hearinghelp.com, by mail: Hearing Help Express • PO Box 586 • DeKalb, IL 60115, or by text: 815-286-4005

Try it for yourself for 45 days risk-free



Order the hearing aid that is right for you. You'll receive high-quality hearing aids with excellent sound that you can **try at home for 45-days risk-free**. If for any reason the model you receive is not quite right, we'll help you find a different option or refund 100% of your purchase price.

Unlimited support to help you every step of the way

Our licensed hearing aid professionals are just one call, video conference or email away, to help you at any point in your journey to better hearing.



Choosing your hearing aid is easy

Hearing aids are not as complicated as some would like you to believe. There are patterns to hearing loss and all hearing aids have one thing in common: they make words louder and easier to understand. The good news is Hearing Help Express® offers you several options to help ***you choose the right price, style, and features you need for improved hearing and understanding!***

Levels of hearing loss symptoms

> Mild to Moderate Hearing Loss Symptoms MOST COMMON

1. Women and children are sometimes hard to understand.
2. You often ask people to repeat themselves.
3. Words seem loud enough, just not clear enough.
4. You can hear in a quiet environment but not in noise.

> Moderate to Severe Hearing Loss Symptoms

1. The above is true and also ...
2. Telephone conversations are often difficult.
3. Any background noise makes it difficult to understand.
4. You sometimes can't hear doorbells or the phone ring.

> Severe to Profound Hearing Loss Symptoms

1. The above is true and also ...
2. You cannot use the telephone without amplification.
3. You often do not hear your alarm clock.
4. You can only understand people when they are looking at you and raise their voice.

To make it easy for you to select the perfect hearing aid, each model has been labeled with the appropriate hearing loss level. If you need any help, we encourage you to call and talk to a hearing aid professional: 1-800-221-2099.



Explanation of hearing aid features

Digital – A digital hearing aid is more sophisticated than its analog predecessor. Digital hearing aids convert natural sounds into coded signals, process and refine the signals, then send the improved sound into the ear. They can be programmed to emphasize sounds in certain frequencies and block out unwanted noise more effectively, allowing for improved hearing and understanding.

Feedback Management – Technology that monitors the hearing aid for feedback (whistling). When detected, the whistling will be reduced or eliminated.

Noise Reduction – Analyzes signals to determine if they contain unwanted noise. If unwanted noise is detected, the hearing aid will reduce the amplification of this unwanted noise, making background noise less noticeable.

Directional Microphones – Designed to enhance sound that is coming from in front of you and reduces amplification of the sounds coming from other directions. This technology improves speech understanding.

Telephone Setting – Also known as a T-Coil or Telecoil. This feature works by turning off the microphone on your hearing aid and only picking up the sounds coming through the phone or loop systems.

The telephone setting is an optional feature that can be added at no additional charge to many of our Apollo and Ranger Series hearing aids as well as our Felix-8750 hearing aid model.

Voice Prompts – This technology may “tell” you when your battery is low or needs to be replaced, what volume level you’re on or what listening setting you are using.



Questions? Call toll-free: 1-800-221-2099
Monday - Friday • 7:30 am - 5:00 pm Central Time

EarMate® Series



Features
Features explained on pg. 7

	EarMate-4220	EarMate-4330	EarMate-4210
Hearing Loss Level	Mild to Moderate	Mild to Moderate	Mild
Hearing Assessment Score	6-25	6-29	6-17
100% Digital	✓	✓	✓
American-Made Technology	✓	✓	✓
Feedback Management	✓	✓	✓
Noise Reduction	✓	✓	✓
Telephone Setting (T-Coil)	X	X	X
Directional Microphones	X	X	X
Voice Prompts	X	X	✓
Volume Control Type	Fingertip	Fingertip	ACOUS-TAP
Battery Size	10	10	10

The EarMate Series features three discreet ready-to-wear in-the-ear hearing aid options. Perfect for those with a mild to moderate hearing loss whether you're new to hearing aids or an experienced user.

- Discreet in-the-ear hearing aids
- Easy to use and maintain
- Perfect for daily use or just when you need a little extra help

**1-800-221-2099
ext. 30601**



EarMate®-4220

- Budget-friendly
- Discreetly fits in the ear
- Simple to use
- Fingertip volume control

Our budget-friendly, EarMate-4220, is a discreet in-the-ear hearing aid designed for mild to moderate hearing loss. This easy-to-use hearing aid is perfect for every day use or just when you need a little extra help. Features a fingertip volume dial and micro-handle for easy inserting and removing.

Average Battery Life: 82 hrs

Color Choice Beige

\$399 each
or \$99.75 per aid for 4 months

HEARING ASSESSMENT SCORE: 6-25
for mild to moderate hearing loss

Customer Reviews

“Comfortable. I can hear the kitchen clock - haven't heard it in years. Family and friends can use a lower voice when talking to me.”

Dennis – Sioux Falls, SD

“Affordable price. I have been wearing your hearing aids for about 15 years. Happy to be able to communicate. Thanks again for your great service.”

William Weiman – PA

“I like the way it fits and how well I can hear with it. Church services are much better. I wear it all the time.”

Harry England – TN



EarMate®-4330

- > Small and easy to use
- > Fingertip volume dial
- > Multiple listening settings

The EarMate-4330 is our most advanced in-the-ear hearing aid. Designed with 100% digital technology, this ready-to-wear hearing aid delivers crisp, clear sound quality for true listening enjoyment. Includes ACOUS-TAP technology that allows you to easily switch between multiple listening settings with a simple tap over your ear. Customize the fit of the EarMate-4330 with soft, flexible eartips.



Average Battery Life: 75 hrs

Color Choice Beige

\$599 each
or \$149.75 per aid for 4 months

HEARING ASSESSMENT SCORE: 6-29
for mild to moderate hearing loss

EarMate®-4210

- > Super discreet
- > Select from multiple volume levels using ACOUS-TAP technology

Designed for those with a busy lifestyle, slip in the EarMate-4210 when you need a boost or wear it all day for steady, comfortable listening. Switch between multiple volume



levels using ACOUS-TAP technology, simply tap on your ear to change the volume. This nearly invisible hearing aid is so small and delivers a crystal clear hearing experience, you might forget you're wearing it!



ACTUAL SIZE



Not recommended for those with dexterity issues.

Average Battery Life: 100 hrs

Color Choice Black

\$699 each
or \$174.75 per aid for 4 months

HEARING ASSESSMENT SCORE: 6-17
for mild hearing loss



Apollo® Series

Features

Features explained on pg. 7

	Apollo-6205	Apollo-E600	Apollo-6670
Hearing Loss Level	Mild to Moderate	Mild to Severe	Mild to Severe
Hearing Assessment Score	6-29	18-41	18-41
100% Digital	✓	✓	✓
American-Made Technology	✓	✓	✓
Feedback Management	✓	✓	✓
Noise Reduction	✓	✓	✓
Telephone Setting (T-Coil)	Optional	X	Optional
Directional Microphones	✓	X	✓
Voice Prompts	X	X	✓
Volume Control Type	Rocker	Wheel	Rocker
Battery Size	312	RECHARGEABLE	312

The Apollo Series features three slim-tube hearing aids. These small, lightweight hearing aids are designed to hide discreetly behind the ear. A very thin tube routes the sound from the hearing aid into the ear. The slim-tubing connects to a soft eartip that sets in the ear canal, giving you the most natural sound with total comfort.

Slim-tube provides extra comfort for those that wear glasses.



**1-800-221-2099
ext. 30610**

Apollo®-6205



- > Budget-friendly
- > Rocker style volume control
- > Multiple listening settings
- > Optional telephone setting

The Apollo-6205 is our budget-friendly slim-tube behind-the-ear hearing aid. The Apollo-6205 is a great starter hearing aid for those with a mild to moderate hearing loss and is also an ideal option to have on hand as a spare. Includes an easy to use rocker style volume control and a push button control that allows you to select from up to four listening settings.

Average Battery Life: 200 hrs Color Choice Beige

\$299 each
or \$74.75 per aid for 4 months

HEARING ASSESSMENT SCORE: 6-29
for mild to moderate hearing loss

Apollo®-E600 *RECHARGEABLE!*

- > Convenient and easy to use
- > fingertip volume control
- > 3 program settings
- > Full charge gives you up to 20 hours of hearing time

Make improved hearing convenient with the rechargeable Apollo-E600 hearing aid. This FDA-Listed hearing aid includes premium digital electronics so you can communicate, and participate more fully in daily activities. For personalized hearing results, select from three program settings and dial in your preferred volume level with the fingertip volume wheel.



Rechargeable: Charging cord included Color Choice Beige

\$699 each
or \$174.75 per aid for 4 months

HEARING ASSESSMENT SCORE: 18-41
for mild to severe hearing loss

Apollo®-6670 ***Best technology for the price!***

- > Advanced digital technology
- > Up to 4 listening settings including an optional telephone setting
- > 10 volume levels



We've taken a customer favorite and given it a boost in technology. The Apollo-6670 processes sound faster to give you more accurate hearing results. This high-end hearing aid fits a wide range of hearing loss and is packed with 100% digital technology that ensures conversations remain clear, even in noisy environments. You'll get improved hearing in a wide range of listening situations encountered in everyday life. **The Apollo-6670 is your best value in terms of technology and price!**

Average Battery Life: 202 hrs Color Choice Bronze or Metallic Black

**\$649 each
or \$162.25 per aid for 4 months**

**HEARING ASSESSMENT SCORE: 18-41
for mild to severe hearing loss**

Customer Reviews

“It works like a charm. Finally I can understand what is said on my TV. Best investment I ever made.”

Carl T – Stockton, CA

“Small, lightweight and easy to use. I have never had a hearing aid that worked this well. Clarity is very impressive. Thank you Hearing Help Express.”

Rick S –Pierceton, IN

“I can hear so much better and people don't have to repeat themselves. This hearing aid is wonderful.”

Edna Elliott –GA

Felix™-8750 ACCU-FIT




So much more than just a hearing aid!


The Felix-8750 Accu-Fit features our very best technology, loaded with personalized programming to give you a one-of-a-kind hearing experience. Also includes Hearing Help Express's Teleaudiology services, providing you with several added benefits via the HHE Hearing app.

The Felix-8750 fits a wide range of hearing loss and is available in the discreet slim-tube style as well as the standard fit behind-the-ear for those that need more power. Designed with premium digital technology including directional microphones that help you focus on the conversation, not the noise around you. The Felix-8750 provides you with the best listening performance in noisy situations.

The Felix-8750 can be completely personalized for your specific hearing loss. Discuss your hearing with one of our licensed professionals or send us a copy of your most recent hearing test. Our professionals will program the Felix-8750 to meet your specific hearing needs.



You are given full control of your Felix hearing aids. Easily adjust the volume and program settings from your smart phone. Plus on-demand support! Simply tap the support button to connect to a licensed hearing aid professional. (iOS and Android.)



Get ready for a one-of-a-kind hearing experience

- ✓ Premium digital technology
- ✓ Personalized for your specific hearing needs
- ✓ Improved hearing in noise
- ✓ Control with your smartphone



Felix-8750 ACCU-FIT

Hearing Loss Level	Mild - Severe
Hearing Assessment Score	6-41
100% Digital	✓
American-Made Technology	✓
Feedback Management	✓
Noise Reduction	✓
Telephone Setting (T-Coil)	Optional
Directional Microphones	✓
Voice Prompts	✓
Companion microphone	Optional
Wireless connectivity	✓
On-demand support	✓
Battery Size	312

Average Battery Life: varies

Color Choice Bronze Metallic Black

Hear the words, not the noise!

The Felix-8750 can be paired wirelessly to the *optional* companion microphone, putting you closer to the sounds you want to hear without all the background noise. Place it next to your TV speaker, clip it to the shirt of a friend or loved one, the Companion Microphone picks up the sound it is closest to and sends it directly to the Felix hearing aid(s).

- **Enjoy one-on-one conversations** without missing a word
- **Hear what's being said** in group conversations, even in noisy environments
- **Hear every word** of your favorite TV or radio program
- **Improved listening** from a distance



Optional Companion Microphone

**Felix-8750 – \$749 each
or \$187.25 per aid for 4 months**

Optional Companion Microphone.....\$100

FREE Companion Microphone when you purchase a pair of Felix-8750 hearing aids!

Ranger® Series



Features Features explained on pg. 7	Ranger-8720	Ranger-7800	Ranger-7910
Hearing Loss Level	Moderate - Severe	Severe - Profound	Severe - Profound
Hearing Assessment Score	30-53	54-60	54-60
100% Digital	✓	✓	✓
American-Made Technology	✓	✓	✓
Feedback Management	✓	✓	✓
Noise Reduction	✓	✓	✓
Telephone Setting (T-Coil)	Optional	✓	✓
Directional Microphones	✓	X	X
Voice Prompts	✓	X	X
Volume Control Type	Push-button	Wheel	Wheel
Battery Size	312	13	675

The Ranger® Series features three standard fit behind-the-ear models. The Ranger Series hearing aids allow for more power with less chance of feedback because sensitive electronic components rest behind the ear. This also makes them less susceptible to earwax damage, requiring less frequent repairs. *Custom fitted earmolds can be added to the Ranger Series for the best fit possible. Call for details.*

Slips behind
the ear. Can
comfortably be
worn with
glasses.



**1-800-221-2099
ext. 30610**

Ranger®-8720

- > More power, more options, lower price
- > Control from your smartphone
- > 3 listening settings and 10 volume levels
- > Optional telephone setting

The Ranger-8720, designed with advanced technology, gives you improved hearing in almost any listening environment. Multiple listening settings and digital programming give you crisp, clear sound.

You have full control with the Ranger-8720. Simply connect the Ranger-8720 to the HHE Hearing app and you

can discreetly make changes to the volume level or listening settings with a few simple taps on your smart phone. The HHE Hearing app is available for both iPhone and Android users.



The Ranger-8720 can be paired wirelessly to the *optional* companion microphone. The companion microphone picks up the sounds it is closest to and sends it directly to the Ranger hearing aid(s). Improves listening from a distance, one-on-one conversations, and when you're in a noisy listening environment.



With the Ranger-8720, you get more power and more options for a great low price.

Average Battery Life: 100 hrs Color Choice Beige

\$599 each
or \$149.75 per aid for 4 months

HEARING ASSESSMENT SCORE: 30-53
for moderate to severe hearing loss

You can purchase the optional Companion Microphone for \$100.

Ranger®-7800



- > Compact and powerful
- > 4 listening settings including telephone setting
- > Fingertip volume control

The Ranger-7800 is one of our smallest high powered hearing aids so if you need a lot of power packed into a compact hearing aid then the Ranger-7800 is right for you. This small but powerful hearing aid features technology that improves overall sound quality and makes speech easier to understand. Includes telephone setting to make phone conversations easier to understand.

Average Battery Life: 195 hrs Color Choice Beige

DUE TO THE HIGH POWER OF THE RANGER-7800, CALL BEFORE ORDERING

**\$899 each
or \$224.75 per aid for 4 months**

**HEARING ASSESSMENT SCORE: 54-60
for severe to profound hearing loss**

Ranger®-7910

- > High-powered for those who need it
- > Long battery life
- > 4 program settings including noise reduction and telephone settings



The Ranger-7910 will help almost anyone with a severe to profound hearing loss. This high-powered hearing aid features 4 listening settings, including a noise setting for improved hearing in noisy environments and feedback cancellation to reduce unwanted squeal. Enjoy the telephone setting that makes your phone conversations easy to hear.

Average Battery Life: 250 hrs Color Choice Beige

DUE TO THE HIGH POWER OF THE RANGER-7910, CALL BEFORE ORDERING

**\$999 each
or \$249.75 per aid for 4 months**

**HEARING ASSESSMENT SCORE: 54-60
for severe to profound hearing loss**

Protect your hearing investment

Accident Protection Plan

Hearing Help Express hearing aids are manufactured to last many years with good care. However, hearing aids are fragile electronic devices that can break with normal wear and tear and more frequently with accidents.

A hearing aid placed near the sink may get wet, your dog may think it's a chew toy, you could drop it, step on it and crack the shell... You'd be surprised how often these things happen.

Avoid costly hearing aid repairs and protect your hearing aid with the Accident Protection Plan. You'll have peace of mind knowing your hearing aids are protected from expensive repairs and common accidents. Coverage includes damage due to earwax & moisture build-up. Loss and theft are not covered.

Accident Protection Plan Pricing

Pricing is per hearing aid*

	<u>1 Year</u>	<u>2 Years</u>
EarMate-4220	\$90	\$158
EarMate-4330	\$110	\$193
EarMate-4210	\$ 95	\$166
Apollo-6205	\$ 90	\$158
Apollo-E600	\$135	\$236
Apollo-6670	\$115	\$201
Felix-8750	\$130	\$228
Ranger-8720	\$125	\$219
Ranger-7800	\$210	\$368
Ranger-7910	\$240	\$420

- > Protection from water damage, droppage, even if you step on it!
- > We'll repair any damage or send a replacement hearing aid at no charge!
- > Simple quick process. We will diagnose and repair or replace within just a few days.

With the Protection Plan, loaner hearing aids are available for only \$15 each (regularly priced at \$40 each).

*Protection plan prices are only valid when purchasing the plan at the same time as the hearing aid.

Cleaning Care Plan – \$149⁹⁹ per hearing aid for one year of coverage

Today's hearing aids are filled with sensitive electronics and over time, earwax and moisture can affect those electronics, weakening the performance of your hearing aids. Choose the Cleaning Care Plan to have your aids cleaned and serviced by our team of in-house technicians to get the most out of your hearing aids.

- **Extend the Life of Your Aids -**

Four professional cleanings per year. Simply call when you're ready for service. Customers who have their hearing aids cleaned by our technicians replace their hearing aids less frequently than those who do not.

- **Our Process -** Up to four times per year we'll send you loaner hearing aids with prepaid postage to send your aids to us for service. Once service is complete we will return your hearing aids along with prepaid postage for you to mail the loaner aids back to us.

Premium Hearing Aid Batteries

Long-Lasting Power • Guaranteed Fresh • Great Low Prices



ONLY
69¢
PER BATTERY

Eco-Gold MAX

REGULAR TABS
Budget friendly favorite

36 Eco-Gold MAX batteries

\$24.84 (69¢ EACH)



ONLY
71¢
PER BATTERY

powerone))

LONGER TABS
Easier to handle

36 Power One batteries

\$25.56 (71¢ EACH)



ONLY
75¢
PER BATTERY

RAYOVAC®

REGULAR TABS
Trusted brand since 1930

40 RAYOVAC batteries

\$29.99 (75¢ EACH)

Need help finding your size? All manufacturers use the same color code system!

Battery Size:
Tab color or color on packaging:

10 Yellow	312 Brown	13 Orange	675 Blue
----------------------------	----------------------------	----------------------------	---------------------------

Helpful Accessories

Premium Battery Tester

Put your mind at ease... test your hearing aid battery and you'll know for sure if it has enough power to run your hearing aid. Simply place your battery in the appropriate size slot, squeeze the blue bar and color coding indicates the remaining power.



Activair Battery Tester.....\$27.99

Travel Battery Tester

This pocket battery tester is small enough to go with you anywhere. Digital tester takes the guesswork out of battery testing. Includes key chain and built-in battery holder compartment.

Key Chain Battery Tester.....\$13.99



Battery Magnet Stick

Hearing aid batteries are small and hard to handle. With the Battery Magnet Stick you can easily insert and remove batteries.

Battery Magnet Stick\$5.95



Battery Care Kit

Worry-free hearing aid battery use.

Battery Care Kit includes:

- Digital Keychain Battery Tester
- Battery Caddie (holds 2 batteries)
- Battery Magnet Stick

All items come in a resealable, clear plastic travel bag.

Battery Care Kit\$21.95



Moisture can attack the sensitive electronics in hearing aids and needlessly shorten their life. Prevent damage by storing your hearing aid in a hearing aid dehumidifier.



BEST VALUE DryAid Kit

Our budget friendly hearing aid dehumidifier is great for every day use and it's rechargeable.

BEST VALUE DryAid Kit \$16.95

Mini Dri-Aid Kit

The rechargeable Mini Dri-Aid Kit is smaller and more compact for easy packing when traveling. Mini Dri-Aid Kit ... \$18.99



Super Dri-Aid Kit

Larger size makes this hearing aid dehumidifier more powerful. The Super Dri-Aid Kit includes pellets that can be recharged and used indefinitely.

Super Dri-Aid Kit \$35.99

PerfectDry Lux – The PerfectDry Lux dehumidifier provides 30 minute rapid drying time. Also provides UV-C disinfection which destroy 99% of microorganisms which are often responsible for itchiness, allergies, and other irritations of the ear canal. Includes USB power supply.

PerfectDry Lux \$90.00



Make daily cleaning easy with premium cleaning tools.



In-the-Ear Kit



Behind-the-Ear Kit

Both kits include:

- Larger, easy to hold cleaning brush
- Handy magnet that makes it easier to pick up batteries
- Individually wrapped hearing aid wipe

Behind-the-Ear kit also includes:

- Wax loop tool to help remove grime
- Specialty cleaning pin for slim-tubing

Specialty In-the-Ear Cleaning Kit\$12.99

Specialty Behind-the-Ear Cleaning Kit\$13.99

Hearing aid clean wipes

Clean and disinfect your hearing aids with this handy pack of Tech-Care Clean Wipes. Ultra fabric wipes are pre-moistened with a hospital grade disinfectant. Alcohol-free and safe to use on all hearing products.

(30 wipes per container)

Tech-Care Clean Wipes\$9.95



Disinfectant towelettes

Start your hearing aid maintenance routine by regularly cleaning and disinfecting your hearing aids with these handy individually wrapped pre-moistened towelettes.



(40 individually-wrapped towelettes)

Hearing Aid Wipes\$11.95

NANOCLEAN

NanoClean Advanced Technology features a stiff threader and a soft, spongy, nylon brush that safely and effectively sweeps away debris and absorbs moisture from hearing aid tubing. “Thread-and-pull” cleaning action gets the job done easily! (20 per container)



NANOCLEAN\$10.99



Remove unwanted earwax

Earwax sometimes builds up in ears. Don't use cotton swabs to clean your ears. Safely remove earwax with this handy kit.

Includes: liquid solution and bulb syringe.

Earwax Removal Kit.....\$10.99

Itchy ears? ProEar by MiraCell helps relieve itchy, irritated ears! It has mild antibacterial, anti-fungal and anti-yeast properties. ProEar is a safe, all-natural solution. Use it without worry to maintain healthy, itch-free ears. (*1/2 ounce bottle*)

ProEar by MiraCell.....\$14.99



Vibrating Alarm Clock

Wake up easily and on time with the Sonic Shaker™ vibrating alarm clock. Powerful vibrating motor will wake the heaviest sleeper. Or choose the extra loud audible alarm and snooze feature. The Sonic Shaker also works great as a travel alarm clock.

Sonic Shaker™

Vibrating Alarm Clock...\$65.99



Order Risk-free!

Try any Hearing Help Express accessory for 30 days risk-free. If you're not 100% satisfied, return all materials within 30 days from the date of delivery, and we'll refund 100% of your purchase price.

Accessory prices and availability are subject to change without notice

See our full accessory product line at HearingHelp.com

ORDERING MADE EASY



an IntriCon company

PO Box 586 • DeKalb, IL 60115-0586

FOR FASTEST SERVICE CALL TOLL-FREE
1-800-221-2099

ORDER ONLINE
HearingHelp.com

ORDER BY MAIL: complete this form and mail it in the envelope provided

PROMO
CODE
30610

Customer information

Dr./ Mr./ Mrs./ Ms. _____

Mailing Address _____

City/State/Zip _____

Phone (required) _____ *Include me in emails with offers and updates.*

Ordered Items

Item Description	Quantity	Price	Total
Hearing Aid Model:	Which Ear? Right / Left / Both		
Felix-8750 single aid order:	<input type="checkbox"/> Companion Microphone	\$100	
<input type="checkbox"/> I am purchasing a pair of Felix-8750's, include my FREE Companion Microphone			
Color choice for Apollo-6670 or Felix-8750 only:	<input type="checkbox"/> Bronze <input type="checkbox"/> Metallic Black	\$0	
<input type="checkbox"/> Add telephone setting (Apollo-6205, Apollo-6670, Felix-8750, Ranger-8720)		\$0	
<input type="checkbox"/> PROTECT MY INVESTMENT – add the one or two year protection plan to each hearing aid being ordered. Pricing on page 19: Aid 1: \$ _____ Aid 2: \$ _____			
Shipping option (shipping fees are non-refundable):	Subtotal _____		
FREE Standard Shipping (Delivery: 7-10 business days).....\$ FREE	Shipping _____		
Priority Shipping (Delivery: 4-5 business days).....\$ 5.95			
USPS Express Mail (Delivery: 1-3 business days).....\$ 24.00	Total _____		

Payment information – 100% Money-Back Guarantee!

Total payment enclosed: \$ _____ Check Credit / Debit Card

Monthly payments: I agree to make 4 interest-free monthly payments for each hearing aid starting today.

First payment enclosed: \$ _____ Check Credit / Debit Card

I understand 4 equal monthly payments will be charged to my card beginning with shipment.

To discuss alternative payment plans, call 1-800-221-2099 (M-F • 8 am - 5 pm CT).

Credit Card Information: Visa MasterCard Discover AMEX

Credit Card # _____ Exp Date: _____

Signature: _____

YOUR MONEY-BACK GUARANTEE FROM HEARING HELP EXPRESS®

Your 2-way assurance of complete satisfaction

45-Day, Risk-Free Trial Money-Back Guarantee

Try any Hearing Help Express® hearing aid for 45 days risk-free. If you are not 100% satisfied, simply return all materials within 45 days from the date of delivery, and we will refund any money paid towards the purchase price. You are responsible for the safe return of all materials in their original condition. Ship insured mail.

One Year Limited Manufacturer's Warranty

Your Hearing Help Express hearing aids are warranted for one year against defects in manufacturer's materials and workmanship. For service, you must send your hearing aids to Hearing Help Express for this warranty to apply. Always mail hearing aids in a sturdy box and ship insured mail.

Important Information for Prospective Hearing Aid Users – Your 45-day home trial will be honored for 60 days starting from the date we ship your order. (The extra 15 days ensures that your home trial will not be shortened by shipping delays.) If hearing aids are accepted for return after 60 days they will be subject to a monthly rental fee of 10% of the regular purchase price of the hearing aid. Hearing aids returned damaged may not be accepted and will be subject to a damage fee. All business is transacted in Illinois and governed by Illinois, FDA, and FTC regulations and the laws of interstate commerce. Monthly payments available, subject to approval. Payments are interest-free when made on schedule. Delinquent accounts are charged 1.4% per month (APR 18.2%).

Call or write if you want to read the User Instructional Brochure for a hearing aid before ordering. You should see a doctor before ordering hearing aids if you have:

- * *History of active drainage from the ear within the previous 90 days*
- * *A history of sudden or rapidly progressive hearing loss within the previous 90 days*
- * *Unilateral hearing loss of sudden or recent onset within the previous 90 days*
- * *Acute or chronic dizziness*
- * *Visible evidence of significant wax accumulation or a foreign body in your ear canal*
- * *Visible congenital or traumatic deformity of the ear*
- * *Pain or discomfort in the ear*
- * *Bone conduction thresholds that exceed air conduction thresholds by 15 decibels or more*

Commonly Asked Questions

How can you sell hearing aids legally? The FDA regulates the hearing aid industry. Historically, the FDA required a waiver to be signed by the consumer before purchasing a hearing aid from a direct mail or online hearing aid company. In addition, Licensed Hearing Aid Dispensers are required to be on staff. In 2016, the FDA dropped the waiver requirement allowing adults (over 18 years old) to access hearing aids without being seen by a hearing care professional.

Do I need a hearing test? A hearing test is not required but can be helpful. You can send it to Hearing Help Express for specific recommendations by one of our licensed hearing aid dispensers.

Do I need to see my doctor before ordering? Federal law allows anyone 18 years and older to obtain a hearing aid without consulting a medical doctor. We suggest you discuss your hearing loss with your doctor to be sure nothing else is causing your hearing difficulties.

Satisfaction Survey Results

Dr. Sergei Kochkin, the world's leading industry expert, performed an independent survey of customers who shopped for hearing aids through audiologists or other local professionals, and customers who bought directly from Hearing Help Express.

His research shows that customers who bought their hearing aids directly from Hearing Help Express were more satisfied and felt they received a higher overall value and more benefits than those buying from a traditional retail professional!



VALUE

Hearing Help Express

79%

Audiologists & Retail

65%

OVERALL SATISFACTION

Hearing Help Express

81%

Audiologists & Retail

78%



**Better hearing has never been
so affordable or convenient.**



PROVEN QUALITY



ACCURATE HEARING RESULTS



BEST VALUE HEARING

1-800-221-2099 • HearingHelp.com
1712 Sycamore Rd, PO Box 586 • DeKalb, IL 60115



HH.364-210601

© Hearing Help Express®, Inc